

For Research Use Only

ERp57/ERp60 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 66423-3-PBS



Basic Information

Catalog Number:

66423-3-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG8741

GenBank Accession Number:

BC014433

GeneID (NCBI):

2923

UNIPROT ID:

P30101

Full Name:

protein disulfide isomerase family A,
member 3

Calculated MW:

505 aa, 57 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2D10C11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50755-1, ERp57/ERp60/PDIA3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66423-2-PBS. Detection antibody: 66423-3-PBS. Standard:Ag8741. Range: 0.098-100 ng/mL.



Cytometric bead array standard curve of MP50755-2, ERp57/ERp60/PDIA3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66423-4-PBS. Detection antibody: 66423-3-PBS. Standard:Ag8741. Range: 0.195-100 ng/mL.