

CSN2 Monoclonal antibody, PBS Only

Catalog Number: 66435-1-PBS

Basic Information

Catalog Number: 66435-1-PBS	GenBank Accession Number: BC012629	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 9318	CloneNo.: 2H10E6
Source: Mouse	UNIPROT ID: P61201	
Isotype: IgG1	Full Name: COP9 constitutive photomorphogenic homolog subunit 2 (Arabidopsis)	
Immunogen Catalog Number: AG1415	Calculated MW: 52 kDa	
	Observed MW: 52 kDa	

Applications

Tested Applications:
WB, Indirect ELISA, IHC, IF

Species Specificity:
human, mouse, rat

Background Information


COP9 is an essential component of the COP9 signalosome complex (CSN), a complex involved in various cellular and developmental processes. The CSN complex is an essential regulator of the ubiquitin (Ubl) conjugation pathway by mediating the deneddylation of the cullin subunits of SCF-type E3 ligase complexes, leading to decrease the Ubl ligase activity of SCF-type complexes such as SCF, CSA or DDB2. The complex is also involved in phosphorylation of p53/TP53, c-jun/JUN, I κ B/ κ B, ITPK1 and IRF8/ICBP, possibly via its association with CK2 and PKD kinases. CSN-dependent phosphorylation of TP53 and JUN promotes and protects degradation by the Ubl system, respectively.

Storage


Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only


Selected Validation Data




Immunofluorescent analysis of (4% PFA) fixed human ovary tumor tissue using CSN2 antibody (66435-1-Ig, Clone: 2H10E6) at dilution of 1:400 and CoraLite@488-Conjugated AffiniPure Goat Anti-Mouse IgG(H+L). This data was developed using the same antibody clone with 66435-1-PBS in a different storage buffer formulation.



NIH/3T3 cells were subjected to SDS PAGE followed by western blot with 66435-1-Ig (CSN2 antibody) at dilution of 1:8000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 66435-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human ovary tumor tissue slide using 66435-1-Ig (CSN2 antibody) at dilution of 1:500 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 66435-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human ovary tumor tissue slide using 66435-1-Ig (COPS2 Antibody) at dilution of 1:500 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 66435-1-PBS in a different storage buffer formulation.