

For Research Use Only

CP110 Monoclonal antibody, PBS Only



Catalog Number: 66448-1-PBS

Basic Information

Catalog Number: 66448-1-PBS	GenBank Accession Number: BC036654	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 9738	CloneNo.: 1A3C12
Source: Mouse	UNIPROT ID: O43303	
Isotype: IgG2a	Full Name: CP110 protein	
Immunogen Catalog Number: AG3489	Calculated MW: 991 aa, 111 kDa	
	Observed MW: 110 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human

Background Information

CP110, also named as CCP110 and KIAA0419, is an 110kDa protein. This gene is different to the gene CEP110 (geneID:11064; CNTRL). CP110 is a centriolar protein that positively regulates centriole duplication while restricting centriole elongation and ciliogenesis. And it acts as a key negative regulator of ciliogenesis in collaboration with CEP97 by capping the mother centriole thereby preventing cilia formation.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 66448-1-Ig (CP110 Antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 66448-1-PBS in a different storage buffer formulation.