

For Research Use Only

IRF4 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 66451-2-PBS



Basic Information

Catalog Number:

66451-2-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG25203

GenBank Accession Number:

BC015752

GeneID (NCBI):

3662

ENSEMBL Gene ID:

ENSG00000137265

UNIPROT ID:

Q15306

Full Name:

IRF 4

Calculated MW:

52 kDa

Purification Method:

Protein G purification

CloneNo.:

2B2C4

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP50012-1, IRF4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66451-2-PBS. Detection antibody: HRP-conjugated 66451-1-PBS. Standard: Ag25203. Range: 0.02-20 ng/mL



Sandwich ELISA standard curve of MP50012-2, IRF4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66451-2-PBS. Detection antibody: HRP-conjugated 66451-3-PBS. Standard: Ag25203. Range: 0.02-20 ng/mL