

For Research Use Only

# DCX Monoclonal antibody, PBS Only (Capture)

Catalog Number: 66462-2-PBS



## Basic Information

Catalog Number:

66462-2-PBS

Concentration:

1mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG4945

GenBank Accession Number:

BC027925

GeneID (NCBI):

1641

UNIPROT ID:

O43602

Full Name:

doublecortin

Calculated MW:

402 aa, 45 kDa

Purification Method:

Protein G purification

CloneNo.:

1F11B11

## Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Sandwich ELISA standard curve of MP50026-1, DCX  
Monoclonal Matched Antibody Pair, PBS Only.  
Capture antibody: 66462-2-PBS. Detection  
antibody: HRP-conjugated 66462-3-PBS. Standard:  
Ag4945. Range: 0.156-20 ng/mL