

For Research Use Only

LCP2 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 66465-3-PBS



Basic Information

Catalog Number:
66465-3-PBS

Source:
Mouse

Isotype:
IgG1

Immunogen Catalog Number:
AG3432

GenBank Accession Number:
BC016618

GeneID (NCBI):
3937

UNIPROT ID:
Q13094

Full Name:
lymphocyte cytosolic protein 2 (SH2 domain containing leukocyte protein of 76kDa)

Calculated MW:
533 aa, 60 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:
2H3A7

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

