

For Research Use Only

# VAMP3 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 66488-2-PBS



## Basic Information

Catalog Number:

66488-2-PBS

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG1156

GenBank Accession Number:

BC005941

GeneID (NCBI):

9341

UNIPROT ID:

Q15836

Full Name:

vesicle-associated membrane protein  
3 (cellubrevin)

Calculated MW:

11 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1E12F9

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

# Selected Validation Data



Cytometric bead array sample test of MP51510-1, VAMP3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66488-1-PBS. Detection antibody: 66488-2-PBS.



Cytometric bead array standard curve of MP51510-1, VAMP3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66488-1-PBS. Detection antibody: 66488-2-PBS. Standard: Ag1156. Range: 1.563-200 ng/mL.