

For Research Use Only

# HRPT2; CDC73 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 66490-2-PBS



## Basic Information

Catalog Number:

66490-2-PBS

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG2959

GenBank Accession Number:

BC014351

GeneID (NCBI):

79577

UNIPROT ID:

Q6P1J9

Full Name:

cell division cycle 73, Paf1/RNA  
polymerase II complex component,  
homolog (S. cerevisiae)

Calculated MW:

531 aa, 61 kDa

Purification Method:

Protein A purification

CloneNo.:

3A2C10

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

## Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array sample test of MP51478-1, CDC73 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66490-2-PBS. Detection antibody: 66490-3-PBS.



Cytometric bead array standard curve of MP51478-1, CDC73 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66490-2-PBS. Detection antibody: 66490-3-PBS. Standard: Ag2959. Range: 3.125-400 ng/mL.