

For Research Use Only

DBC1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 66497-2-PBS



Basic Information

Catalog Number:

66497-2-PBS

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG18374

GenBank Accession Number:

BC018269

GeneID (NCBI):

57805

UNIPROT ID:

Q8N163

Full Name:

KIAA1967

Calculated MW:

923 aa, 103 kDa

Purification Method:

Protein G purification

CloneNo.:

2A10A7

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51483-1, CCAR2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66497-2-PBS. Detection antibody: 66497-1-PBS.



Cytometric bead array standard curve of MP51483-1, CCAR2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66497-2-PBS. Detection antibody: 66497-1-PBS. Standard:Ag18374. Range: 1.563-100 ng/mL.