

For Research Use Only

# TAU Monoclonal antibody, PBS Only (Detector)

Catalog Number: 66499-2-PBS



## Basic Information

Catalog Number:

66499-2-PBS

Source:

Mouse

Isotype:

IgG3

Immunogen Catalog Number:

AG21926

GenBank Accession Number:

BC000558

GeneID (NCBI):

4137

UNIPROT ID:

P10636

Full Name:

microtubule-associated protein tau

Calculated MW:

37-46, 79-81 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

1D7C9

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human, mouse, rat

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array sample test of MP51432-1, TAU Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60995-1-PBS. Detection antibody: 66499-2-PBS.