

For Research Use Only

WWTR1/TAZ Monoclonal antibody, PBS Only (Capture)

Catalog Number: 66500-2-PBS



Basic Information

Catalog Number: 66500-2-PBS	GenBank Accession Number: BC014052	Purification Method: Protein G purification
Source: Mouse	GeneID (NCBI): 25937	CloneNo.: 2D9A9
Isotype: IgG1	UNIPROT ID: Q9GZV5	
Immunogen Catalog Number: AG13330	Full Name: WW domain containing transcription regulator 1	
	Calculated MW: 44 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51507-1, TAZ Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66500-2-PBS. Detection antibody: 66500-3-PBS.



Cytometric bead array standard curve of MP51507-1, TAZ Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66500-2-PBS. Detection antibody: 66500-3-PBS. Standard: Ag13330. Range: 1.563-200 ng/mL.