

For Research Use Only

NRF2, NFE2L2 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 66504-4-PBS

Basic Information

Catalog Number: 66504-4-PBS	GenBank Accession Number: BC011558	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 4780	CloneNo.: 1E10H2
Source: Mouse	UNIPROT ID: Q16236	
Isotype: IgG1	Full Name: nuclear factor (erythroid-derived 2)- like 2	
Immunogen Catalog Number: AG9469	Calculated MW: 605 aa, 68 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50316-3, NRF2, NFE2L2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66504-4-PBS. Detection antibody: 66504-5-PBS. Standard:Ag9469. Range: 0.391-25 ng/mL.