

For Research Use Only

Recoverin Monoclonal antibody, PBS Only (Capture)

Catalog Number: 66521-2-PBS



Basic Information

Catalog Number:

66521-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG0119

GenBank Accession Number:

BC001720

GeneID (NCBI):

5957

UNIPROT ID:

P35243

Full Name:

recoverin

Calculated MW:

24 kDa

Purification Method:

Protein G purification

CloneNo.:

1E6A8

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51037-1, Recoverin Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66521-2-PBS. Detection antibody: 66521-3-PBS. Standard:Ag0119. Range: 0.391-100 ng/mL



Cytometric bead array standard curve of MP51037-2, Recoverin Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66521-2-PBS. Detection antibody: 66521-4-PBS. Standard:Ag0119. Range: 1.563-100 ng/mL