

For Research Use Only

# PCCA Monoclonal antibody, PBS Only (Capture)



Catalog Number: 66538-2-PBS

## Basic Information

Catalog Number:

66538-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG17271

GenBank Accession Number:

BC000140

GeneID (NCBI):

5095

UNIPROT ID:

P05165

Full Name:

propionyl Coenzyme A carboxylase,  
alpha polypeptide

Calculated MW:

80 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1D6D5

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP50299-1, PCCA Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66538-2-PBS. Detection antibody: 66538-3-PBS. Standard:Ag17271. Range: 3.125-200 ng/mL.