

For Research Use Only

PCCA Monoclonal antibody, PBS Only (Detector)

Catalog Number: 66538-5-PBS



Basic Information

Catalog Number:

66538-5-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG17271

GenBank Accession Number:

BC000140

GeneID (NCBI):

5095

UNIPROT ID:

P05165

Full Name:

propionyl Coenzyme A carboxylase,
alpha polypeptide

Calculated MW:

80 kDa

Purification Method:

Protein G purification

CloneNo.:

1A12H12

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50299-2, PCCA Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66538-4-PBS. Detection antibody: 66538-5-PBS. Standard: Ag17271. Range: 6.25-200 ng/mL



Biolayer interferometry (BLI) kinetic assay of 66538-5-Ig against Human PCCA was performed. The affinity constant is 5.47 nM.