

For Research Use Only

RYR1 Monoclonal antibody, PBS Only

Catalog Number: 66539-1-PBS



Basic Information

Catalog Number:

66539-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG25120

GenBank Accession Number:

NM_000540

GeneID (NCBI):

6261

Full Name:

ryanodine receptor 1 (skeletal)

Calculated MW:

565 kDa

Purification Method:

Protein A purification

CloneNo.:

1B4E12

Applications

Tested Applications:

Indirect ELISA, IHC

Species Specificity:

Human, Mouse

Background Information

Ryanodine receptor isoform-1 (RyR1) is a major calcium channel in skeletal muscle important for excitation-contraction coupling. The transmembrane region of RyR1 is comprised of two domains including the pseudo voltage sensor domain and the pore-forming domain (PMID: 27855725). Dominant RyR1 mutations are the leading cause of MH. Studies on the Ca²⁺ conducting properties of RyR1 channels containing Malignant hyperthermia mutations expressed in myotubes.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded mouse skeletal muscle tissue slide using 66539-1-Ig (RYR1 antibody) at dilution of 1:200 (under 10x lens. Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 66539-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded mouse skeletal muscle tissue slide using 66539-1-Ig (RYR1 antibody) at dilution of 1:200 (under 40x lens. Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 66539-1-PBS in a different storage buffer formulation.