

For Research Use Only

Calponin Monoclonal antibody, PBS Only (Capture)



Catalog Number: 66540-2-PBS

Basic Information

Catalog Number: 66540-2-PBS	GenBank Accession Number: BC022015	Purification Method: Protein G purification
Size: 1mg/ml	GeneID (NCBI): 1264	CloneNo.: 1B3C9
Source: Mouse	UNIPROT ID: P51911	
Isotype: IgG1	Full Name: calponin 1, basic, smooth muscle	
Immunogen Catalog Number: AG21384	Calculated MW: 33 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP50017-1, Calponin Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66540-2-PBS. Detection antibody: HRP-conjugated 66540-1-PBS. Standard: Ag21384. Range: 0.313-20 ng/mL.