

For Research Use Only

TOP2A Monoclonal antibody, PBS Only (Capture)

Catalog Number: 66541-2-PBS



Basic Information

Catalog Number:	GenBank Accession Number:	Purification Method:
66541-2-PBS	BC140791	Protein G Magarose purification
Size:	GeneID (NCBI):	CloneNo.:
1 mg/ml	7153	1A2A2
Source:	UNIPROT ID:	
Mouse	P11388	
Isotype:	Full Name:	
IgG1	topoisomerase (DNA) II alpha 170kDa	
Immunogen Catalog Number:	Calculated MW:	
AG19551	1531 aa, 174 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51073-1, TOP2A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66541-2-PBS. Detection antibody: 66541-3-PBS. Standard:Ag19551. Range: 0.195-100 ng/mL



Cytometric bead array standard curve of MP51073-4, TOP2A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66541-2-PBS. Detection antibody: 66541-7-PBS. Standard:Ag19551. Range: 0.195-100 ng/mL