

For Research Use Only

# TOP2A Monoclonal antibody, PBS Only (Capture)

Catalog Number: 66541-5-PBS



## Basic Information

Catalog Number:

66541-5-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG19551

GenBank Accession Number:

BC140791

GeneID (NCBI):

7153

UNIPROT ID:

P11388

Full Name:

topoisomerase (DNA) II alpha 170kDa

Calculated MW:

1531 aa, 174 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2G6C10

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP51073-3, TOP2A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66541-5-PBS. Detection antibody: 66541-6-PBS. Standard:Ag19551. Range: 0.195-100 ng/mL.