

For Research Use Only

# HuR Monoclonal antibody, PBS Only (Capture)

Catalog Number: 66549-2-PBS



## Basic Information

Catalog Number:

66549-2-PBS

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG2249

GenBank Accession Number:

BC003376

GeneID (NCBI):

1994

UNIPROT ID:

Q15717

Full Name:

ELAV (embryonic lethal, abnormal vision, Drosophila)-like 1 (Hu antigen R)

Calculated MW:

326 aa, 36 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1G5D9

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array sample test of MP51567-1, HuR Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66549-2-PBS. Detection antibody: 66549-3-PBS.



Cytometric bead array standard curve of MP51567-1, HuR Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66549-2-PBS. Detection antibody: 66549-3-PBS. Standard: Ag2249. Range: 1.563-200 ng/mL