

For Research Use Only

Ki-67 Monoclonal antibody

Catalog Number: 66555-1-Ig **3 Publications**



Basic Information

Catalog Number: 66555-1-Ig	GenBank Accession Number: NM_002417	Purification Method: Protein G purification
Source: Mouse	GeneID (NCBI): 4288	CloneNo.: 4G8A6
Isotype: IgG1	UNIPROT ID: P46013	
Immunogen Catalog Number: AG26266	Full Name: antigen identified by monoclonal antibody Ki-67	
	Calculated MW: 359 kDa	

Applications

Tested Applications:
ELISA

Cited Applications:
IHC, IF

Species Specificity:
Human

Cited Species:
human, rat

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Xiaozhi Li	36316517	J Cancer Res Clin Oncol	IHC
Chaoqiong Ding	41747254	Cancer Discov	IF
Shaojun Fang	40368885	Light Sci Appl	IHC

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data