

For Research Use Only

Ki-67 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 66555-4-PBS



Basic Information

Catalog Number:

66555-4-PBS

Concentration:

1000 µg/ml

Source:

Mouse

Isotype:

IgM

Immunogen Catalog Number:

AG26266

GenBank Accession Number:

NM_002417

GeneID (NCBI):

4288

UNIPROT ID:

P46013

Full Name:

antigen identified by monoclonal antibody Ki-67

Calculated MW:

359 kDa

Purification Method:

Caprylic acid/ammonium sulfate precipitation

CloneNo.:

1A6D2

Applications

Tested Applications:

IHC, Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50237-2, KI67 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66555-4-PBS. Detection antibody: 66555-3-PBS. Standard:Ag26266. Range: 0.391-50 ng/mL