

For Research Use Only

# Ki-67 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 66555-5-PBS



## Basic Information

Catalog Number:

66555-5-PBS

Concentration:

1000 µg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG26266

GenBank Accession Number:

NM\_002417

GeneID (NCBI):

4288

UNIPROT ID:

P46013

Full Name:

antigen identified by monoclonal antibody Ki-67

Calculated MW:

359 kDa

Purification Method:

Protein A purification

CloneNo.:

1D8E4

## Applications

Tested Applications:

IHC, Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP50237-3, KI67 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66555-5-PBS. Detection antibody: 66555-6-PBS. Standard:Ag26266. Range: 0.098-50 ng/mL