

For Research Use Only

OMA1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 66572-1-PBS



Basic Information

Catalog Number: 66572-1-PBS	GenBank Accession Number: BC012915	Purification Method: Protein G purification
Size: 1mg/ml	GeneID (NCBI): 115209	CloneNo.: 2E11D7
Source: Mouse	UNIPROT ID: Q96E52	
Isotype: IgG1	Full Name: OMA1 homolog, zinc metallopeptidase (S. cerevisiae)	
Immunogen Catalog Number: AG11201	Calculated MW: 524 aa, 60 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50571-1, OMA1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66572-1-PBS. Detection antibody: 66572-2-PBS. Standard:Ag11201. Range: 0.781-100 ng/mL



Cytometric bead array standard curve of MP50571-3, OMA1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66572-1-PBS. Detection antibody: 67449-1-PBS. Standard:Ag11201. Range: 3.125-100 ng/mL