

For Research Use Only

# GRP78/BIP Monoclonal antibody, PBS Only (Detector)

Catalog Number: 66574-3-PBS



## Basic Information

Catalog Number:

66574-3-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG2188

GenBank Accession Number:

BC020235

GeneID (NCBI):

3309

UNIPROT ID:

P11021

Full Name:

heat shock 70kDa protein 5 (glucose-regulated protein, 78kDa)

Calculated MW:

78 kDa

Purification Method:

Protein G purification

CloneNo.:

1E5E6

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

# Selected Validation Data



Sample test of MP50855-1, GRP78,BIP Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66574-2-PBS. Detection antibody: 66574-3-PBS. Cell lysate(30  $\mu$ g/well): HepG2, HT-1080, Hela.



Cytometric bead array standard curve of MP50855-1, GRP78,BIP Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66574-2-PBS. Detection antibody: 66574-3-PBS. Standard:Ag2188. Range: 1.563-25 ng/mL.