

For Research Use Only

CHCHD5 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 66584-2-PBS



Basic Information

Catalog Number:

66584-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG22534

GenBank Accession Number:

BC004498

GeneID (NCBI):

84269

UNIPROT ID:

Q9BSY4

Full Name:

coiled-coil-helix-coiled-coil-helix domain containing 5

Calculated MW:

110 aa, 12 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

3F3B4

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50619-1, CHCHD5 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66584-2-PBS. Detection antibody: 66584-3-PBS. Standard:Ag22534. Range: 0.098-50 ng/mL



Cytometric bead array standard curve of MP50619-2, CHCHD5 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66584-2-PBS. Detection antibody: 66584-4-PBS. Standard:Ag22534. Range: 0.098-3.125 ng/mL