

For Research Use Only

# TOM70 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 66593-2-PBS



## Basic Information

Catalog Number:

66593-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG6027

GenBank Accession Number:

BC052994

GeneID (NCBI):

9868

UNIPROT ID:

O94826

Full Name:

translocase of outer mitochondrial membrane 70 homolog A (S. cerevisiae)

Calculated MW:

67 kDa

Purification Method:

Protein A purification

CloneNo.:

2B6B10

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP50612-1, TOM70 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66593-2-PBS. Detection antibody: 66593-3-PBS. Standard:Ag6027. Range: 0.781-100 ng/mL



Cytometric bead array standard curve of MP50612-2, TOM70 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66593-2-PBS. Detection antibody: 66593-4-PBS. Standard:Ag6027. Range: 1.563-100 ng/mL