

For Research Use Only

SAMD3 Monoclonal antibody, PBS Only



Catalog Number: 66599-1-PBS

Basic Information

Catalog Number:

66599-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG11319

GenBank Accession Number:

BC029851

GeneID (NCBI):

154075

UNIPROT ID:

Q8N6K7

Full Name:

sterile alpha motif domain containing

3

Calculated MW:

221aa, 25 kDa; 520aa, 61 kDa

Observed MW:

52 kDa

Purification Method:

Protein G purification

CloneNo.:

1A11F2

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

Human, mouse, pig, rat

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 66599-1-Ig (SAMD3 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 66599-1-PBS in a different storage buffer formulation.