

For Research Use Only

PPP3CB Monoclonal antibody, PBS Only (Capture)

Catalog Number: 66614-2-PBS



Basic Information

Catalog Number: 66614-2-PBS	GenBank Accession Number: BC028049	Purification Method: Protein G purification
Source: Mouse	GeneID (NCBI): 5532	CloneNo.: 1E5A2
Isotype: IgG1	UNIPROT ID: P16298	
Immunogen Catalog Number: AG26987	Full Name: protein phosphatase 3 (formerly 2B), catalytic subunit, beta isoform	
	Calculated MW: 525 aa, 59 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51504-1, PPP3CB Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66614-2-PBS. Detection antibody: 66614-3-PBS.



Cytometric bead array standard curve of MP51504-1, PPP3CB Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66614-2-PBS. Detection antibody: 66614-3-PBS. Standard: Ag26987. Range: 0.781-100 ng/mL.