

For Research Use Only

# PPP3CB Monoclonal antibody, PBS Only (Detector)

Catalog Number: 66614-3-PBS



## Basic Information

Catalog Number:

66614-3-PBS

Source:

Mouse

Isotype:

IgG3

Immunogen Catalog Number:

AG26987

GenBank Accession Number:

BC028049

GeneID (NCBI):

5532

UNIPROT ID:

P16298

Full Name:

protein phosphatase 3 (formerly 2B),  
catalytic subunit, beta isoform

Calculated MW:

525 aa, 59 kDa

Purification Method:

Protein G purification

CloneNo.:

1C10H7

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array sample test of MP51504-1, PPP3CB Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66614-2-PBS. Detection antibody: 66614-3-PBS.



Cytometric bead array standard curve of MP51504-1, PPP3CB Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66614-2-PBS. Detection antibody: 66614-3-PBS. Standard: Ag26987. Range: 0.781-100 ng/mL.