

For Research Use Only

RABEPK Monoclonal antibody, PBS Only (Detector)

Catalog Number: 66622-2-PBS



Basic Information

Catalog Number:

66622-2-PBS

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG7796

GenBank Accession Number:

BC065725

GeneID (NCBI):

10244

UNIPROT ID:

Q7Z6M1

Full Name:

Rab9 effector protein with kelch motifs

Calculated MW:

41 kDa

Purification Method:

Protein G purification

CloneNo.:

1F11A3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51484-1, RABEPK Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66622-1-PBS. Detection antibody: 66622-2-PBS.



Cytometric bead array standard curve of MP51484-1, RABEPK Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66622-1-PBS. Detection antibody: 66622-2-PBS. Standard:Ag2593. Range: 1.563-400 ng/mL