

For Research Use Only

# GRIK2 Monoclonal antibody, PBS Only



Catalog Number: 66631-2-PBS

## Basic Information

Catalog Number:

66631-2-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG25308

GenBank Accession Number:

BC037954

GeneID (NCBI):

2898

UNIPROT ID:

Q13002

Full Name:

glutamate receptor, ionotropic,  
kainate 2

Calculated MW:

908 aa, 103 kDa

Observed MW:

103-115 kDa

Purification Method:

Protein G purification

CloneNo.:

1E5G6

## Applications

Tested Applications:

Indirect ELISA, WB

Species Specificity:

pig, rat, mouse, human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 66631-2-Ig (GRIK2 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 66631-2-PBS in a different storage buffer formulation.