

For Research Use Only

ETV5 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 66657-1-PBS



Basic Information

Catalog Number: 66657-1-PBS	GenBank Accession Number: BC007333	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 2119	CloneNo.: 2A4A4
Source: Mouse	UNIPROT ID: P41161	
Isotype: IgG2a	Full Name: ets variant 5	
Immunogen Catalog Number: AG3708	Calculated MW: 510 aa, 58 kDa	
	Observed MW: 58-60 kDa	

Applications

Tested Applications:
WB, Sandwich ELISA, Indirect ELISA

Species Specificity:
Human

Background Information

ETS variant gene 5(ETV5) belongs to the ETS oncogene family that shares a conserved peptide, ETS domain, which mediates sequence-specific DNA binding. ETV5 is a transcription factor and is required for spermatogonial stem cell self-renewal. ETV5 can be negatively regulated by COP1, a tumor suppressor. Also it has a role in branching morphogenesis in the developing kidney

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 66657-1-Ig (ETV5 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 66657-1-PBS in a different storage buffer formulation.



Sandwich ELISA standard curve ofMP50024-1, ETV5 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66657-1-PBS. Detection antibody: HRP-conjugated 66657-2-PBS. Standard: Ag3708. Range: 0.313-20 ng/m.