

For Research Use Only

ETV5 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 66657-2-PBS

Basic Information

Catalog Number:

66657-2-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG3708

GenBank Accession Number:

BC007333

GeneID (NCBI):

2119

UNIPROT ID:

P41161

Full Name:

ets variant 5

Calculated MW:

510 aa, 58 kDa

Purification Method:

Protein G purification

CloneNo.:

4E3C2

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve ofMP50024-1, ETV5
Monoclonal Matched Antibody Pair, PBS Only.
Capture antibody: 66657-1-PBS. Detection
antibody: HRP-conjugated 66657-2-PBS. Standard:
Ag3708. Range: 0.313-20 ng/m.