

For Research Use Only

STK4 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 66663-3-PBS



Basic Information

Catalog Number: 66663-3-PBS	GenBank Accession Number: BC093768	Purification Method: Protein G purification
Source: Mouse	GeneID (NCBI): 6789	CloneNo.: 3D1H11
Isotype: IgG1	UNIPROT ID: Q13043	
Immunogen Catalog Number: AG17738	Full Name: serine/threonine kinase 4	
	Calculated MW: 487 aa, 56 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Sample test of MP51388-1, STK4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66663-2-PBS. Detection antibody: 66663-3-PBS.



Cytometric bead array standard curve of MP51388-1, STK4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66663-2-PBS. Detection antibody: 66663-3-PBS. Standard: Ag17738. Range: 0.098-100 ng/mL



Bi-layer interferometry (BLI) kinetic assays of 66663-3-Ig against Human STK4 were performed. The affinity constant is 245 nM.