

For Research Use Only

# MLKL Monoclonal antibody, PBS Only (Capture)

Catalog Number: 66675-2-PBS



## Basic Information

Catalog Number: 66675-2-PBS	GenBank Accession Number: BC028141	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 197259	CloneNo.: 2F8E7
Source: Mouse	UNIPROT ID: Q8NB16	
Isotype: IgG1	Full Name: mixed lineage kinase domain-like	
Immunogen Catalog Number: AG15188	Calculated MW: 471 aa, 54 kDa	

## Applications

Tested Applications:  
Cytometric bead array, Indirect ELISA

Species Specificity:  
human

## Background Information

### Storage

Storage:  
Store at -80°C.  
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP51035-1, MLKL Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66675-2-PBS. Detection antibody: 66675-3-PBS. Standard:Ag15188. Range: 6.25-100 ng/mL



Cytometric bead array standard curve of MP51035-2, MLKL Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66675-2-PBS. Detection antibody: 66675-4-PBS. Standard:Ag15188. Range: 3.125-100 ng/mL