

For Research Use Only

HMG20A Monoclonal antibody, PBS Only (Detector)

Catalog Number: 66676-3-PBS



Basic Information

Catalog Number:

66676-3-PBS

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG2718

GenBank Accession Number:

BC021959

GeneID (NCBI):

10363

UNIPROT ID:

Q9NP66

Full Name:

high-mobility group 20A

Calculated MW:

347 aa, 40 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

1C4D3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51552-1, HMG20A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66676-2-PBS. Detection antibody: 66676-3-PBS.



Cytometric bead array standard curve of MP51552-1, HMG20A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66676-2-PBS. Detection antibody: 66676-3-PBS. Standard:Ag2718. Range: 3.125-400 ng/mL.