

For Research Use Only

# ERC1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 66681-3-PBS



## Basic Information

Catalog Number: 66681-3-PBS	GenBank Accession Number: BC132782	Purification Method: Protein G purification
Source: Mouse	GeneID (NCBI): 23085	CloneNo.: 2A4B6
Isotype: IgG1	UNIPROT ID: Q8IUD2	
Immunogen Catalog Number: AG17665	Full Name: ELKS/RAB6-interacting/CAST family member 1	
	Calculated MW: 1116 aa, 128 kDa	

## Applications

Tested Applications:  
Cytometric bead array, Indirect ELISA

Species Specificity:  
human

## Background Information

### Storage

Storage:  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:  
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

# Selected Validation Data



Cytometric bead array sample test of MP51549-2, ERC1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66681-3-PBS. Detection antibody: 66681-4-PBS.



Cytometric bead array standard curve of MP51549-2, ERC1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66681-3-PBS. Detection antibody: 66681-4-PBS. Standard: Ag17665. Range: 0.781-100 ng/mL.