

For Research Use Only

ABCD3 Monoclonal antibody, PBS Only



Catalog Number: 66697-1-PBS

Basic Information

Catalog Number:

66697-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG3

Immunogen Catalog Number:

AG24760

GenBank Accession Number:

BC068509

GeneID (NCBI):

5825

UNIPROT ID:

P28288

Full Name:

ATP-binding cassette, sub-family D (ALD), member 3

Observed MW:

70 kDa

Purification Method:

Protein A purification

CloneNo.:

1G7G5

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 66697-1-Ig (ABCD3 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 66697-1-PBS in a different storage buffer formulation.