

For Research Use Only

TPT1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 66713-3-PBS



Basic Information

Catalog Number: 66713-3-PBS	GenBank Accession Number: BC003352	Purification Method: Protein G purification
Source: Mouse	GeneID (NCBI): 7178	CloneNo.: 2A7F7
Isotype: IgG1	UNIPROT ID: P13693	
Immunogen Catalog Number: AG27651	Full Name: tumor protein, translationally-controlled 1	
	Calculated MW: 20 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage


Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3


For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51548-1, TPT1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66713-2-PBS. Detection antibody: 66713-3-PBS.



Cytometric bead array standard curve of MP51548-1, TPT1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66713-2-PBS. Detection antibody: 66713-3-PBS. Standard: Ag27651. Range: 3.125-400 ng/mL.