

For Research Use Only

ARNT/HIF-1 beta Monoclonal antibody, PBS Only (Capture)

Catalog Number: 66732-2-PBS



Basic Information

Catalog Number:

66732-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG5507

GenBank Accession Number:

BC060838

GeneID (NCBI):

405

UNIPROT ID:

P27540

Full Name:

aryl hydrocarbon receptor nuclear
translocator

Calculated MW:

416 aa, 46 kDa, 87 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

4E1E8

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50351-1, ARNT,HIF1B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66732-2-PBS. Detection antibody: 66732-3-PBS. Standard:Ag5507. Range: 0.195-6.25 ng/mL.



Cytometric bead array standard curve of MP50351-2, ARNT,HIF1B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66732-2-PBS. Detection antibody: 66732-4-PBS. Standard:Ag5507. Range: 0.781-25 ng/mL.