

For Research Use Only

TDP-43 Monoclonal antibody

Catalog Number: 66734-2-Ig



Basic Information

Catalog Number: 66734-2-Ig	GenBank Accession Number: BC001487	Purification Method: Protein G purification
Size: 1000 µg/ml	GeneID (NCBI): 23435	CloneNo.: 3C5H2
Source: Mouse	UNIPROT ID: Q13148	
Isotype: IgG1	Full Name: TAR DNA binding protein	
Immunogen Catalog Number: AG1231	Calculated MW: 43 kDa	

Applications

Tested Applications:
ELISA

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

