

For Research Use Only

Androgen receptor Monoclonal antibody, PBS Only (Capture)

Catalog Number: 66747-6-PBS



Basic Information

Catalog Number:

66747-6-PBS

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG17291

GenBank Accession Number:

BC132975

GeneID (NCBI):

367

UNIPROT ID:

P10275

Full Name:

androgen receptor

Calculated MW:

914 aa, 99 kDa

Purification Method:

Protein A purification

CloneNo.:

2C1F9

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51133-3, Androgen receptor Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66747-6-PBS. Detection antibody: 66747-5-PBS. Standard:Ag17291. Range: 0.391-200 ng/mL.



Biolayer interferometry (BLI) kinetic assay of 66747-6-Ig against Human Androgen receptor was performed. The affinity constant is below 1 pM.