

For Research Use Only

KIAA0284 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 66749-2-PBS

Basic Information

Catalog Number:

66749-2-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG26852

GenBank Accession Number:

BC112928

GeneID (NCBI):

283638

UNIPROT ID:

Q9Y4F5

Full Name:

KIAA0284

Calculated MW:

1589 aa, 172 kDa

Purification Method:

Protein G purification

CloneNo.:

1G4E3

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP50021-1, KIAA0284 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66749-2-PBS. Detection antibody: HRP-conjugate 66749-1-PBS. Standard: Ag26852. Range: 0.313-20 ng/mL.