

For Research Use Only

# Collagen Type I Monoclonal antibody, PBS Only (Capture)

Catalog Number: 66761-2-PBS



## Basic Information

Catalog Number:

66761-2-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG6281

GenBank Accession Number:

BC054498

GeneID (NCBI):

1278

UNIPROT ID:

P08123

Full Name:

collagen, type I, alpha 2

Calculated MW:

1366 aa, 130 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1H9F4

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP51074-1, COL1A2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66761-2-PBS. Detection antibody: 66761-3-PBS. Standard:Ag6281. Range: 0.098-100 ng/mL.