

For Research Use Only

Apolipoprotein A II/APOA2 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 66773-2-PBS



Basic Information

Catalog Number:

66773-2-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG9863

GenBank Accession Number:

BC005282

GeneID (NCBI):

336

UNIPROT ID:

P02652

Full Name:

apolipoprotein A-II

Calculated MW:

100 aa, 11 kDa

Purification Method:

Protein A purification

CloneNo.:

2C1C3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50949-1, APOA2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66773-2-PBS. Detection antibody: 66773-3-PBS. Standard:Ag9863. Range: 0.098-100 ng/mL.