

For Research Use Only

# STX16 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 66775-2-PBS



## Basic Information

Catalog Number:

66775-2-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG27926

GenBank Accession Number:

BC019042

GeneID (NCBI):

8675

UNIPROT ID:

O14662

Full Name:

syntaxin 16

Calculated MW:

37 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

3H5B4

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Sample test of MP51331-1, STX16 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66775-1-PBS. Detection antibody: 66775-2-PBS.



Cytometric bead array standard curve of MP51331-1, STX16 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66775-1-PBS. Detection antibody: 66775-2-PBS. Standard:Ag27926. Range: 1.563-200 ng/mL