

For Research Use Only

CD52 Monoclonal antibody

Catalog Number: 66784-1-Ig



Basic Information

Catalog Number:

66784-1-Ig

Size:

1000 µg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG17788

GenBank Accession Number:

BC000644

GeneID (NCBI):

1043

UNIPROT ID:

P31358

Full Name:

CD52 molecule

Calculated MW:

6 kDa

Observed MW:

20 kDa

Purification Method:

Protein A purification

CloneNo.:

1B10E8

Recommended Dilutions:

WB 1:1000-1:4000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human

Positive Controls:

WB : human spleen tissue,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



human spleen tissue were subjected to SDS PAGE followed by western blot with 66784-1-Ig (CD52 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.