

For Research Use Only

RGS3 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 66790-3-PBS



Basic Information

Catalog Number:

66790-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG3177

GenBank Accession Number:

BC018072

GeneID (NCBI):

5998

UNIPROT ID:

P49796

Full Name:

regulator of G-protein signaling 3

Calculated MW:

1198 aa, 132 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1F1B4

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50850-1, RGS3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66790-2-PBS. Detection antibody: 66790-3-PBS. Standard:Ag3177. Range: 1.563-200 ng/mL.