

For Research Use Only

AQP1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 66805-3-PBS



Basic Information

Catalog Number:

66805-3-PBS

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG14093

GenBank Accession Number:

BC022486

GeneID (NCBI):

358

UNIPROT ID:

P29972

Full Name:

aquaporin 1 (Colton blood group)

Calculated MW:

269 aa, 29 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

2B5F7

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51521-1, AQP1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66805-2-PBS. Detection antibody: 66805-3-PBS.



Cytometric bead array standard curve of MP51521-1, AQP1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66805-2-PBS. Detection antibody: 66805-3-PBS. Standard:Ag14093. Range: 1.563-200 ng/mL.